




# Transport Vials

30ml




**FOR TRANSPORT OF  
FAECAL SPECIMENS**

FOR INTESTINAL  
PARASITE EXAMINATION



## PARASITOLOGY


Single Use In Vitro Diagnostic Device



# Transport Vials

## 30ml

### PRODUCT ORDERING



In the laboratory intestinal parasites are confirmed by the identification of cysts, protozoan trophozoites, larvae and helminth eggs. The priorities of these clinical laboratories do not always permit the immediate examination of a fresh sample.

The prompt collection and transportation of these samples cannot always be guaranteed.

Furthermore, freezing, refrigeration and / or incubation of samples cannot ensure the full recovery of parasites at all stages of identification.

The ReIngenuity Transportation Vial will preserve the intestinal parasites in the faecal material until such time it can be examined by a qualified parasitologist.

#### Product Description

#### Pack Size

30ml Transport Vial (Clean)	40
30ml Transport Vial with 15ml Formalin	40
30ml Transport Vial with 15ml SAF	40
30ml Transport Vial with 15ml Universal Fixative	40



Reingenuity Limited  
Brentwood Units, Lincoln Road,  
Market Rasen,  
LN8 3SF  
info@reingenuity.com

